



CERABIEN FULL KIT

CERABIEN is a porcelain specially developed to make all ceramic crowns (bridges), using Alumina frameworks. It does not cause any allergies and has the advantage to have no black lines in the cervical area. The thermal expansion of CERABIEN is totally different from that of porcelain fused to metal. Therefore, mixing or using with porcelain fused to metal is not recommended.

Coefficient of thermal expansion
(50–500°C 10⁻⁶ K⁻¹)

EX-3	Cerabien	CZR
12.4	6.8	9.1

SHADE BASE	10g each	18 colors
BODY		18 colors
ENAMEL		3 colors
TRANSLUCENT		4 colors
CERVICAL		4 colors
MODIFIER		14 colors
ADD-ON		2 colors
COLOR COMBINATION TABLE		1 piece
TECHNICAL INSTRUCTIONS		1 copy

CERABIEN MARGIN KIT

MARGIN PORCELAIN is developed exclusively for the combination application with Alumina framework. In addition to ordinary usage, it can be used to correct a chip in the Alumina framework. The use of MAGIC FORMER as forming liquid is highly recommended.

BASIC COLOR	10g each	14 colors
SPECIAL COLOR		3 colors
DILUTION POWDER (MDL)		1 color
RETOUCHING POWDER (MRP)		1 color

CERABIEN OPACIOUS BODY KIT

OPACIOUS BODY can be used to achieve delicate colors and to control translucency.

BASIC COLOR	10g each	16 colors
SPECIAL COLOR		3 colors

CERABIEN LUSTER KIT



LUSTER PORCELAIN reproduces natural enamel in both surface features and color.

- LUSTER PORCELAIN reproduces the fine surface structure and luster of natural teeth.
- A unique combination of fine surface particles produces a selective reflection of light that results in the same opalescence seen in natural teeth.
- LUSTER PORCELAIN has transparent, bright and vivid colors, therefore darkening at the incisal edge or at the occlusal surface will not occur.
- Noritake has thoroughly studied color changes in natural teeth caused by aging. LUSTER PORCELAIN features a complete line of colors consistent with these changes.
- **“Creamy White” is also available as individual item.**

BASIC COLOR	10g each	4 colors
SPECIAL COLOR		3 colors

CERABIEN INTERNAL STAIN KIT



Using INTERNAL STAIN, it is possible to replicate many of the most intricate three-dimensional patterns of shade variation found in natural teeth.

1. WHITE	7. CERVICAL 2	13. A+	3g each
2. INCISAL BLUE 1	8. CERVICAL 3	14. B+	
3. INCISAL BLUE 2	9. EARTH BROWN	15. C+	
4. MAWELON ORANGE 1	10. REDDISH BROWN	16. D+	
5. MAWELON ORANGE 2	11. SALMON PINK		
6. CERVICAL 1	12. RED		
BRIGHT (DILUTION)		3g	1 color
IS LIQUID		10ml	1 bottle
IS COLOR GUIDE			1 piece
TECHNICAL INSTRUCTIONS			1 copy

CERABIEN FORMING LIQUID

Specially formulated for CERABIEN.

MEISTER LIQUID

This liquid is slow-drying forming liquid. Keep porcelain moist, which makes applying Shade Base or other powders easily.

IS LIQUID

Does not change the color of INTERNAL STAIN after being mixed, which allows the confirmation of colors before baking.

ES LIQUID

Keeps EXTERNAL STAIN moist, which makes staining work easy.

CERABIEN EXTERNAL STAIN KIT



EXTERNAL STAIN reproduces the great variety of color shades observed on the surface of natural teeth.

1. PURE WHITE	8. ORANGE 1	15. SALMON PINK	3g each
2. GRAY	9. ORANGE 2	16. PINK	
3. BLACK	10. CERVICAL 1	17. RED	
4. BLUE	11. CERVICAL 2	18. A+	
5. GREEN 1	12. CERVICAL 3	19. B+	
6. GREEN 2	13. EARTH BROWN	20. C+	
7. YELLOW	14. REDDISH BROWN	21. D+	
GLAZE		10g	1 color
ES LIQUID		10ml	1 bottle
ES COLOR GUIDE			1 piece
TECHNICAL INSTRUCTIONS			1 copy