



CZR-SS FULL KIT

### CZR FULL KIT / CZR-SS FULL KIT

CZR FULL KIT is an excellent kit composed of basic colors for all ceramic crowns and bridges, using Zirconia frameworks.

CZR-SS FULL KIT has the same shade selection as CZR FULL KIT. Shade Base Stains are contained instead of Shade Base porcelains.

Coefficient of thermal expansion  
(50–500°C 10<sup>-6</sup> K<sup>-1</sup>)

EX-3	Cerabien	CZR
12.4	6.8	9.1

SB (powder) / SS (stain)	10g each (SS 3g each)	18 colors
BODY		18 colors
ENAMEL		3 colors
TRANSLUCENT		4 colors
CERVICAL		4 colors
MODIFIER		14 colors
ADD-ON		2 colors
COLOR COMBINATION TABLE		1 piece
TECHNICAL INSTRUCTIONS		1 copy



### CZR MARGIN KIT

In addition to ordinary usage, MARGIN PORCELAIN can be used to correct a chip in the Zirconia framework. The use of MAGIC FORMER as forming liquid is highly recommended.

BASIC COLOR	10g each	14 colors
SPECIAL COLOR		3 colors
DILUTION POWDER (MDL)		1 color
RETOUCHING POWDER (MRP)		1 color



### CZR OPACIOUS BODY KIT

OPACIOUS BODY can be used to achieve delicate colors and to control translucency.

BASIC COLOR	10g each	16 colors
SPECIAL COLOR		3 colors



### CZR LUSTER KIT

LUSTER PORCELAIN reproduces natural enamel in both surface features and color.

- LUSTER PORCELAIN reproduces the fine surface structure and luster of natural teeth.
- A unique combination of fine surface particles produces a selective reflection of light that results in the same opalescence seen in natural teeth.
- LUSTER PORCELAIN has transparent, bright and vivid colors, therefore darkening at the incisal edge or at the occlusal surface will not occur.
- Noritake has thoroughly studied color changes in natural teeth caused by aging. LUSTER PORCELAIN features a complete line of colors consistent with these changes.

LUSTER COLOR	10g each	8 colors
--------------	----------	----------

### CZR PROFESSIONAL KIT / CZR-SS PROFESSIONAL KIT



CZR-SS PROFESSIONAL KIT

CZR PROFESSIONAL KIT is suitable for expert technicians. It contains a variety of types and colors, including clear cervical porcelain and new internal stain. Now you can approach a natural appearance ever more closely!

CZR-SS PROFESSIONAL KIT has Shade Base Stains instead of Shade Base porcelains.

SB (powder) /SS (stain)		SB White/SS White
ENAMEL	10g each	Silky E1, Silky E2
LUSTER TRANSLUCENT	(SS 3g each)	ELT1, ELT2, ELT3
CLEAR CERVICAL		CCV-1, CCV-2, CCV-3, CCV-4
INTERNAL STAIN	3g	Gray

### CZR ESTHETIC WHITE KIT / CZR-SS ESTHETIC WHITE KIT



CZR-SS ESTHETIC WHITE KIT

CZR ESTHETIC WHITE KIT is for creating whiter shades rather than NW shades. It is not only pure bleaching shades but also natural and beautiful white shades.

CZR-SS ESTHETIC WHITE KIT has a Shade Base Stain (SSB1) instead of Shade Base SBB1.

SB (powder) /SS (stain)		SB White, SBB1/SS White, SSB1
BODY	10g each	EWooB, EWoB, EWB, EWYB
ENAMEL	(SS 3g each)	Silky E1, Silky E2
LUSTER TRANSLUCENT		ELT1, ELT2, ELT3

### CZR INTERNAL STAIN KIT



Using INTERNAL STAIN, it is possible to replicate many of the most intricate three-dimensional patterns of shade variation found in natural teeth.

1. WHITE	7. CERVICAL 2	13. A+	3g each
2. INCISAL BLUE 1	8. CERVICAL 3	14. B+	
3. INCISAL BLUE 2	9. EARTH BROWN	15. C+	
4. MAMELON ORANGE 1	10. REDDISH BROWN	16. D+	
5. MAMELON ORANGE 2	11. SALMON PINK		
6. CERVICAL 1	12. RED		
BRIGHT (DILUTION)	3g	1 color	
IS LIQUID	10ml	1 bottle	
IS COLOR GUIDE		1 piece	
TECHNICAL INSTRUCTIONS		1 copy	

### CZR EXTERNAL STAIN KIT



EXTERNAL STAIN reproduces the great variety of color shades observed on the surface of natural teeth.

1. PURE WHITE	8. ORANGE 1	15. SALMON PINK	3g each
2. GRAY	9. ORANGE 2	16. PINK	
3. BLACK	10. CERVICAL 1	17. RED	
4. BLUE	11. CERVICAL 2	18. A+	
5. GREEN 1	12. CERVICAL 3	19. B+	
6. GREEN 2	13. EARTH BROWN	20. C+	
7. YELLOW	14. REDDISH BROWN	21. D+	
GLAZE	10g	1 color	
ES LIQUID	10ml	1 bottle	
ES COLOR GUIDE		1 piece	
TECHNICAL INSTRUCTIONS		1 copy	

### CZR FORMING LIQUID

Specially formulated for CZR.

100ml

### IS LIQUID

Specially formulated for Internal Stain. It can be used for CZR PRESS Shade Base Stain & CZR PRESS Glaze powder.

10ml

### MEISTER LIQUID

This liquid is slow-drying forming liquid, which gives more time for porcelain application.

100ml

### ES LIQUID

Keeps EXTERNAL STAIN moist, which makes staining work easy.

10ml