



NORITAKE SUPER PORCELAIN

EX-3 PRESS TRIAL KIT L1 / L2

EX-3 PRESS TRIAL KIT L1 & L2 are an introductory kit which includes EX-3 PRESS Opaque, L ingots, EX-3 PRESS LF Enamel, LF Translucent, LF Luster, AD-T, AD-B and MRP Porcelains. TRIAL KIT L1 contains the following 5 shades: NW₀, NW_{0.5}, nA₁, nA₂ & nA₃. And TRIAL KIT L2 contains the following 5 shades: nA₂, nA₃, nA_{3.5}, nB₂ & nC₂. For EX-3 PRESS procedures, Noritake PRESS Investment, Ring / former, Paste Opaque Liquid & ES Liquid are required.

TRIAL KIT L1

PPnA ₁	PPnA ₂	PPnA ₃	PPNW ₀
PPNW _{0.5}	LnA ₁	LnA ₂	LnA ₃
LNW ₀	LNW _{0.5}	LF Silky E ₂	LF E ₂
LF E ₃	LF LT ₁	LF CE	LF Tx
LF AD-T	LF MRP		
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			1 piece

TRIAL KIT L2

PPnA ₂	PPnA ₃	PPnA _{3.5}	PPnB ₂
PPnC ₂	LnA ₂	LnA ₃	LnA _{3.5}
LnB ₂	LnC ₂	LF E ₂	LF E ₃
LF LT ₁	LF CE	LF Tx	LF AD-T
LF AD-B	LF MRP		
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 piece



NORITAKE SUPER PORCELAIN

EX-3 PRESS TRIAL KIT H1 / H2

EX-3 PRESS TRIAL KIT H1 & H2 are an introductory kit which includes EX-3 PRESS Opaque, H ingots, EX-3 PRESS LF External Stain A⁺ & EX-3 PRESS LF Glaze Powder. TRIAL KIT H1 contains the following 5 shades: NW₀, NW_{0.5}, nA₁, nA₂ & nA₃. And TRIAL KIT H2 contains the following 5 shades : nA₂, nA₃, nA_{3.5}, nB₂ & nC₂. The beginner EX-3 PRESS user should also use Noritake PRESS Investment, Ring / former, Paste Opaque Liquid & ES Liquid.

TRIAL KIT H1

HnA ₁	HNW ₀	PPnA ₂	PPNW _{0.5}
HnA ₂	HNW _{0.5}	PPnA ₃	A ⁺
HnA ₃	PPA ₁	PPNW ₀	GLAZE
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces

TRIAL KIT H2

HnA ₂	HnC ₂	PPnB ₂	C ⁺
HnA ₃	PPnA ₂	PPnC ₂	GLAZE
HnA _{3.5}	PPn ₃	A ⁺	
HnB ₂	PPn _{3.5}	B ⁺	
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces

NORITAKE SUPER PORCELAIN

NORITAKE ACT LIQUID



This is Anti-Greening liquid for EX-3 Press material. Sometimes it is necessary to use ACT liquid for any alloy that contains some silver (Ag) to prevent greening.

ACT LIQUID	20ml
------------	------

NORITAKE SUPER PORCELAIN

EX-3 PRESS DISPO PLUNGER

NORITAKE EX-3 PRESS DISPO PLUNGER is disposable plunger set for EX-3 PRESS. 2 kinds of size are available for each 2g & 5g EX-3 PRESS INGOT.



NORITAKE SUPER PORCELAIN **EX-3 PRESS INGOT KIT L**

NORITAKE EX-3 PRESS INGOT KIT L contains low translucency pressable ceramic ingots. It was developed to recreate the margin area precisely by pressing INGOT L. The L ingot is indicated for use in the layering technique. After ingot L is pressed, the crown can be made with EX-3 PRESS LF Enamel, LF Translucent or Luster porcelain, LF Internal stain and LF External stain.

LnA1	LnA2	LnA3	LnA3.5
LnA4	LnB1	LnB2	LnB3
LnB4	LnC1	LnC2	LnC3
LnC4	LnD2	LnD3	LnD4
LNW0	LNW0.5	LNP1.5	LNP2.5
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces



NORITAKE SUPER PORCELAIN **EX-3 PRESS INGOT KIT H**

NORITAKE EX-3 PRESS INGOT-H Kit contains high translucency pressable ceramic ingots. It was developed to recreate both the margin area, as well as the full-contour form of the tooth. The H ingot is indicated for use in the staining method. After pressing the H ingot, the restoration may be characterized using EX-3 PRESS LF External Stain & EX-3 PRESS Glaze Powder. The EX-3 PRESS LF Layering Assortment may also be used after pressing.

HnA1	HnB1	HnC2	HnD4
HnA2	HnB2	HnC3	HNW0
HnA3	HnB3	HnC4	HNW0.5
HnA3.5	HnB4	HnD2	HNP1.5
HnA4	HnC1	HnD3	HNP2.5
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces



NORITAKE SUPER PORCELAIN **EX-3 PRESS PASTE OPAQUE BASIC KIT**

EX-3 PRESS PASTE OPAQUE BASIC KIT is specially designed for EX-3 PRESS INGOT. It can make the basic shade color of metal framework with a thinner coat before pressing ingot.

PPnA1	PPnB1	PPnC2	PPnD4
PPnA2	PPnB2	PPnC3	PPNW0
PPnA3	PPnB3	PPnC4	PPNW0.5
PPnA3.5	PPnB4	PPnD2	PPNP1.5
PPnA4	PPnC1	PPnD3	PPNP2.5
PPEW0	PPEW		
TECHNICAL INSTRUCTIONS			1 copy
PASTE OPAQUE MODIFIER BRUSH			1 piece



NORITAKE SUPER PORCELAIN **EX-3 PRESS PASTE OPAQUE MODIFIER KIT**

NORITAKE EX-3 PRESS PASTE OPAQUE MODIFIER KIT contains 8 color modifiers which may be incrementally incorporated into EX-3 PRESS Basic Paste Opaques to create custom shades. They may also be applied individually for minor color adjustments & subtle shade variances.

1. PP WHITE	5. PP REDDISH BROWN	3g each
2. PP GRAY	6. PP PINK	
3. PP ORANGE	7. PP BLUE	
4. PP EARTH BROWN	8. PP YELLOW	
TECHNICAL INSTRUCTIONS		1 copy
PASTE OPAQUE MODIFIER BRUSH		1 piece



NORITAKE SUPER PORCELAIN **EX-3 PRESS ESTHETIC WHITE KIT**

EX-3 PRESS Esthetic White Kit is for creating whiter shades rather than NW shades, which is not only pure bleaching shades but also natural and beautiful white shades.

PPEW0	PPEW	EW00	EW0
EW	EWY	LF Silky E1	LF Silky E2
LF ELT1	LF ELT2	LF ELT3	
TECHNICAL INSTRUCTIONS			1 copy
EX-3 PRESS DISPO PLUNGER			3 pieces



NORITAKE SUPER PORCELAIN **EX-3 PRESS LF BODY KIT**

EX-3 PRESS LF body can be used for repairing and for additional touch-up on pressable area. This Body Kit is composed of All n Colors, EW, NW and NP shades.

LF nA1B	LF nA2B	LF nA3B	LF nA3.5B
LF nA4B	LF nB1B	LF nB2B	LF nB3B
LF nB4B	LF nC1B	LF nC2B	LF nC3B
LF nC4B	LF nD2B	LF nD3B	LF nD4B
LF EW00B	LF EW0B	LF EWB	LF EWYB
LF NW0B	LF NW0.5B	LF NP1.5B	LF NP2.5B



NORITAKE SUPER PORCELAIN **EX-3 PRESS LF ASSORTMENT KIT**

EX-3 PRESS LF is low fusing porcelain to be built up not only to make an enamel layer but also cervical layer after pressing EX-3 press ingot. This kit contains LF Enamel, Translucent, Luster, Clear Cervical, Modifier, Add-on, MRP porcelains. By using these porcelains with pressable ceramic, you can make more natural crown.

LF E1	LF E2	LF E3	LF Tx	LF T0
LF T1	LF T2	LF LT0	LF LT1	LF T-BLUE
LF LT NATURAL	LF CREAMY ENAMEL	LF SUN BRIGHT	LF INCISAL AUREOLA	LF CREAMY WHITE
LF CCV-1	LF CCV-2	LF CCV-3	LF CCV-4	LF Mamelon 1
LF Mamelon 2	LF Tissue 1	LF Tissue 2	LF Tissue 3	LF Tissue 4
LF AD-T	LF AD-B	LF MRP		
TECHNICAL INSTRUCTIONS			1 copy	

NORITAKE SUPER PORCELAIN **EX-3 PRESS LF INTERNAL STAIN KIT**



EX-3 PRESS LF INTERNAL STAIN KIT is Low fusing stain for EX-3 PRESS LF porcelain. It is necessary to use this stain before EX-3 PRESS LF layering

LF White	LF Incisal Blue 1	LF Incisal Blue 2	LF Mamelon Orange 1
LF Mamelon Orange 2	LF Cervical 1	LF Cervical 2	LF Cervical 3
LF Earth Brown	LF Reddish Brown	LF Salmon Pink	LF Red
LF A+	LF B+	LF C+	LF D+
LF Gray	LF Bright		
TECHNICAL INSTRUCTIONS			1 copy
IS Liquid 10ml			1 bottle
IS Color Guide			1 piece

NORITAKE SUPER PORCELAIN **EX-3 PRESS LF EXTERNAL STAIN KIT**



NORITAKE EX-3 PRESS LF EXTERNAL STAIN KIT is specifically & exclusively formulated for external staining of EX-3 PRESS & EX-3 PRESS stand-alone, pressable restorations. EX-3 LF External Stains enable the ceramist to reproduce the complex surface colors, stains & subtle characterizations visible on the surface of natural teeth.

PURE WHITE	YELLOW	EARTH BROWN	B+	3g each
GRAY	ORANGE 1	REDDISH BROWN	C+	
BLACK	ORANGE 2	SALMON PINK	D+	
BLUE	CERVICAL 1	PINK		
GREEN 1	CERVICAL 2	RED		
GREEN 2	CERVICAL 3	A+		
GLAZE			10g	1 color
ES LIQUID			10ml	1 bottle
ES COLOR GUIDE				1 piece
TECHNICAL INSTRUCTIONS				1 copy