

Super Porcelain Ti-22 DENTAL PORCELAIN for Pure Titanium



Ti-22



Company A

Ti-22 can produce bright colors without being grayish.

NORITAKE SUPER PORCELAIN Ti-22 FULL KIT

NORITAKE SUPER PORCELAIN Ti-22 is for fusing to pure titanium. It has excellent bonding strength because of its Bonding Porcelain that was specially developed for that purpose. It has excellent acid resistance so it is very chemically stable in the mouth. Bright vivid colors are easy to achieve. It does not have the grayish tinge that is common in low fusing porcelains.

OPAQUE	10g each	16 colors
BODY		16 colors
ENAMEL		3 colors
TRANSLUCENT		3 colors
CERVICAL		4 colors
ADD-ON		2 colors
BONDING PORCELAIN		1 color
MODIFIER		12 colors
OPAQUE MODIFIER		8 colors
FORMING LIQUID		10ml each
OPAQUE LIQUID	1 bottle	
BP LIQUID	1 bottle	
COLOR COMBINATION TABLE		1 piece
TECHNICAL INSTRUCTIONS		1 copy



NORITAKE SUPER PORCELAIN Ti-22 MARGIN PORCELAIN KIT

Ti-22 MARGIN PORCELAIN is to be used in combination with Ti-22 PORCELAIN, which is for use with pure titanium. It has low shrinkage to facilitate making margins with excellent fit.

BASIC COLOR	10g each	13 colors
MARGIN MODIFIER		3 colors
POWDER FOR DILUTION		1 color
POWDER FOR RETOUCHING		1 color
COLOR TABLE		1 piece



NORITAKE SUPER PORCELAIN Ti-22 LUSTER PORCELAIN SET

This is for use with NORITAKE SUPER PORCELAIN Ti-22. It has the same features as EX-3 LUSTER PORCELAIN. It is for expert technicians who seek a very natural appearance and complex shades.

1. T BLUE	5. SUN BRIGHT	10g each
2. LUSTER To	6. INCISAL AUREOLA	
3. LUSTER T1	7. CREAMY WHITE	
4. CREAMY ENAMEL		
COLOR TABLE		1 piece
TECHNICAL INSTRUCTIONS		1 copy



NORITAKE SUPER PORCELAIN Ti-22 STAIN KIT

This kit contains external stain for use with NORITAKE SUPER PORCELAIN Ti-22. Because the kit has many colors, it is easy to achieve detailed characterization.

1. PURE WHITE	7. RED	13. CERVICAL 1	3g each
2. GRAY	8. GREEN 1	14. CERVICAL 2	
3. BLACK	9. GREEN 2	15. CERVICAL 3	
4. BLUE	10. GREEN 3	16. EARTH BROWN	
5. YELLOW	11. ORANGE 1	17. REDDISH BROWN	
6. PINK	12. ORANGE 2		
GLAZING POWDER		5g	1 bottle
STAIN LIQUID		10ml	1 bottle
COLOR GUIDE			1 piece

Ti-22 BONDING LIQUID

Best suited for mixing with Ti-22 BONDING PORCELAIN because it shows the ideal compatibility with Titanium.

10ml

Ti-22 OPAQUE LIQUID

Developed for mixing Ti-22 OPAQUE Porcelain in order to burn out at its baking temperature.

10ml 30ml

Ti-22 FORMING LIQUID

Burns out completely at the baking temperature of Ti-22 BODY and ENAMEL porcelain.

10ml 100ml

Ti-22 STAIN LIQUID

Keeps Ti-22 STAIN moist and best suited for mixing with Ti-22 STAIN porcelain.

10ml